

icare
HOME

TA022-088-3.1


EN	INSTRUCTION MANUAL
DE	BEDIENUNGSANLEITUNG
ES	MANUAL DE INSTRUCCIONES
IT	MANUALE DI ISTRUZIONI
SV	ANVÄNDARHANDBOK
FI	KÄYTTÖOHJE
NO	BRUKSANVISNING
PL	INSTRUKCJA OBSŁUGI
DA	BRUGERMANUAL
PT	MANUAL DE INSTRUÇÕES
NL	GEBRUIKSAANWIJZING
FR	MODE D'EMPLOI

1


2

3


4

5

6

7

8

9


10

11

12

13

14


15


16

17


18


19


20


21

22

23

24


25


26

27

28

29


30

31

32

33


34


35


36


37


38

39


40

41

42

43

44

45

ENGLISH

1.	Safety instructions.....	1
2.	Glossary.....	3
3.	Indications for use.....	3
4.	How does the tonometer work	3
5.	Contraindications.....	4
6.	Risks.....	4
7.	Benefits.....	4
8.	Package contents	4
9.	Before you start.....	4
10.	Setting up the tonometer.....	5
	10.1 Installing or changing the battery	5
	10.2 Turning the tonometer on	5
	10.3 Loading the probe.....	5
11.	Using the tonometer.....	5
	11.1 Choosing the measurement mode	5
	11.2 Automatic eye recognition	6
	11.3 Adjusting the measurement position.....	6
	11.4 Taking the measurements	6
12.	Transferring the measurement data	7
13.	Troubleshooting.....	7
14.	Replacing the probe base.....	8
15.	Cleaning and disinfection.....	9
16.	Accessories, detachable parts and other supplies.....	9
17.	Lifetime.....	9
18.	Technical and performance data	9
19.	Symbols	10
20.	Electromagnetic declaration	11
21.	Training.....	13

ENGLISH


Icare® HOME (Model: TA022) INSTRUCTION MANUAL

The information in this document is subject to change without prior notice. In a conflict situation the English version prevails.

This device complies with: Medical Device Directive 93/42/EEC
Canadian Medical Device Regulations
RoHS Directive 2011/65/EU

Copyright © 2020 Icare Finland Oy
Made in Finland

 Icare Finland Oy
Äyritie 22, 01510 Vantaa
Tel. +358 9 8775 1150, Fax +358 9 728 6670
www.icaretonometer.com, info@icarefinland.com

1. SAFETY INSTRUCTIONS

WARNING!

Do not push the tonometer into the eye (the tip of the probe should be 4-8 mm, or 5/32-5/16", from the eye). Never use a probe that has no plastic tip.

WARNING!

Keep the tonometer out of the reach of children and pets, because the probe base, battery compartment cover and probes are so small that a child could swallow them.

WARNING!

The probes are for single-use only, and are packed sterile.

WARNING!

To prevent contamination, do not touch the bare probe, do not use a probe if it touches a non-sterile surface like a table or a floor.

WARNING!

Health care professionals must inform patients not to modify or discontinue their treatment plan without receiving instructions from the health care professional.

WARNING!

For cybersecurity do not connect to the USB port except when uploading patient measurement data. Also the tonometer does not allow you to take any measurements when the USB is connected. All other traditional cybersecurity controls (anti-virus software, malware software, separate network for the device, etc.) do not apply since the device is stand alone, is not networked, and does not contain operating system software.

WARNING!

Do not change the batteries or probe base when the USB cable is connected.

WARNING!

No modification of this equipment is allowed.

WARNING!

Use only the original and certified probes made by the manufacturer. The probes are for single-use (single pair of measurement sequences) only. Use probes taken only from the intact, original packaging. The manufacturer cannot guarantee sterility of the probe once the seal is compromised. Re-sterilization or re-use of the probe could result in incorrect measurement values, in the breakdown of the probe, cross-contamination of bacteria or viruses, and infection of the eye. Re-sterilization or re-use will void all responsibilities and liabilities of the manufacturer concerning the safety and effectiveness of the tonometer.

WARNING!

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

⚠ WARNING!

Never immerse the Icare tonometer in liquid. Do not spray, pour or spill liquid onto the Icare tonometer, its accessories, connectors, switches or openings in the chassis. Remove any liquid appearing on the surface of the tonometer at once.

⚠ WARNING!

The tonometer should only be opened by qualified service personnel. It contains no user-serviceable parts, apart from the batteries and a probe base. The Icare tonometer requires no routine servicing or calibration other than changing the batteries at least every 12 months and the probe base every 6 months. If there is reason to believe servicing of the device is necessary, contact qualified service personnel or your local Icare representative.

⚠ WARNING!

In the event of battery short circuit, the batteries and the device handle could get hot and there is a possible risk of a burn if touched. In such case the device handle could heat up to 60 °C and batteries up to 90 °C and the device should not be used.

⚠ WARNING!

When reading the measurement data in a clinic or hospital environment make sure the tonometer and the PC/mobile device are located a minimum of 1.5 m from the patient.

⚠ WARNING!

Do not use the tonometer in the restricted environments defined in Chapter 18 of this manual.

⚠ WARNING!

Do not connect anything to the USB port but the USB cable supplied with the tonometer.

⚠ WARNING!

For assistance, if needed, in setting up, using or maintaining the tonometer or to report unexpected operation or events, contact the manufacturer or the manufacturer's representative.

⚠ WARNING!

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

⚠ WARNING!

Sources of power frequency magnetic field should be used no closer than 15 cm (6 inches) to any part of Icare HOME (TA022), including cables specified by the manufacturer. Otherwise, degradation of the performance could result.

i PRECAUTION!

- When you have opened the package, check for any external damage or faults, particularly for damage to the case. If you suspect that there is something wrong with the tonometer, contact the dealer who sold the tonometer to you.
- Use the tonometer only for measuring intraocular pressure. Any other use is improper and the manufacturer is not liable for any damage arising from improper use, or for the consequences of such use.
- Never open the casing of the tonometer, except for the battery compartment.
- Never allow the tonometer to get wet.
- Do not use the tonometer near flammable substances, including flammable anesthetic agents.
- Certain microbiological agents (for example, bacteria) can be transmitted from the forehead or cheek support. To prevent this, clean the forehead and cheek support for each new patient with disinfectant. See the chapter 'Cleaning and disinfection'.
- The tonometer conforms to EMC requirements (IEC 60601-1-2), but interference may occur within the tonometer if used near (<1 m) a device causing high intensity electromagnetic emissions, such as a cellular phone. Although the tonometer's own electromagnetic emissions are well below

the levels permitted by the relevant standards, they may cause interference in other, nearby devices, for example sensitive sensors.

- If you do not use the tonometer for a long time, remove the batteries, as they may leak.
- Be sure to dispose of the single-use probes properly (for example, in a mixed waste bin).
- Batteries, packaging materials and probe bases must be disposed of according to local regulations.
- Make sure you use batteries with built-in PTC protection, for example Energizer Lithium Photo 123 3V CR123A.
- Do not cover the eye recognition transmitters or sensor during the measurement, for example with your fingers. Keep your hand, hair etc. and objects such as pillows away from the temple side of your eye, as they produce an infrared reflection that causes an error.
- The tonometer turns off automatically after 3 minutes if you do not use it.
- Update the tonometer's time to your local time. It is done automatically by performing the steps 1 and 2 under section 12. Reading the measurement data.
- Check that the silicone lid covers the USB port during the measurement.
- The measurement method of the Icare HOME tonometer is based on magnetic induction and therefore an external magnetic field in line with the probe may prevent the measurement. In such case the tonometer will continuously ask to repeat the measurement. Situation can be solved either by removing the source of interference from the vicinity of the device or by performing the measurement in different location with no such interference.

2. GLOSSARY

Tonometer = a device for measuring your IOP

mmHg = units of measurement for the pressure in your eye

Position tag = support positions for your eyes are written on it

Probe = tonometer's single-use item that lightly touches your eye

Probe base light = light ring that assists you to position the probe

Shelf life = the time the probe remains sterile in its intact packing

Cornea = the eye's outermost dome shaped clear layer

Forehead/cheek supports = tonometer's adjustable supports

Expected service life = expected life time before replacement

3. INDICATIONS FOR USE

The Icare HOME tonometer is intended for monitoring of intraocular pressure (IOP) of the human eye. It is indicated for use by patients or their caregivers under supervision of an eye care professional.

4. HOW DOES THE TONOMETER WORK

The tonometer measures your eye pressure using a disposable probe that gently contacts your eye during the measurement. You may need a mirror to help you determine the correct measurement position. The recommended frequency of measurements is according to the healthcare professional's instruction or 3-6 times per day.

A complete measurement is a series of six very rapid measurements. It may take from several seconds to a minute to take the measurement. The probe moves to the cornea and back during each of the six rapid measurement. After the series of six measurements, the tonometer calculates your final eye pressure and stores it in the tonometer's memory.

The probe is disposable. You can use the same probe for both eyes if the health care professional instructed you to measure both eyes. After you have taken the measurements for both eyes place the probe back in its container and dispose of it in a mixed waste bin.

The tonometer includes infra-red eye recognition sensors to identify which eye, right or left, you are measuring. It is important not to cover these sensors with your fingers, hand, hair etc., because covering the sensors causes an error. It is also important to keep objects such as pillows away from the temple side of your eye, as they produce an infrared reflection that causes an error.

Do not take the device into use before you have been trained in its use. Do not modify or discontinue your treatment plan without receiving guidance from a health care professional.

The Essential Performance of the Icare HOME (TA022) tonometer is to measure the intraocular pressure (IOP) with specified accuracy (section 18) and transfer the data (section 12).

5. CONTRAINDICATIONS

You should not use the HOME tonometer if you:

1. Have an active ocular infection (for example, pink eye)
2. Had recent trauma to your eye
3. Have disabling arthritis or difficulty handling the Icare tonometer
4. Have a severe difficulty in opening your eyes including abnormal contractions or twitches of the eyelid
5. Have involuntary, rapid and repetitive movement of the eyes

6. RISKS

You may not be eligible for use of the HOME tonometer if you:

1. Have poor uncorrected near vision
2. Have only one working eye
3. Have poor or off-center visual fixation
4. Have poor hearing without an aid and/or communicate using sign language
5. Use contact lens
6. Have dry eyes
7. Have Keratoconus (a corneal disorder)
8. Have a congenitally (from birth) small eye
9. Have an enlarged eye from childhood glaucoma

The safety and effectiveness of the Icare HOME tonometer has not been evaluated for patients with:

1. High corneal astigmatism
2. History of prior invasive glaucoma surgery or corneal surgery including corneal laser surgery (for example, LASIK)
3. Corneal scarring
4. Very thick or very thin corneas
5. Patients for whom there is already difficulty in obtaining clinical IOP measurements (for example, due to eyelid squeezing or tremor)

7. BENEFITS

You can take measurements of your eye pressure at various times during the day and night with the Icare HOME tonometer. These measurements taken outside your office visits may help your doctor better understand your eye pressure, especially if you already have a glaucoma diagnosis.

8. PACKAGE CONTENTS

The package contains:

- Icare HOME tonometer
- 10 sterilized single-use probes
- 2 batteries
- USB memory stick including the instruction manual
- USB cables
- Support position tags
- Warranty card
- Carrying case
- Wrist strap
- Spare probe base

9. BEFORE YOU START

Find the main parts, buttons and indicator lights of the tonometer. The figures are located at the beginning of this document.

FRONT PARTS (Figure 1)

1. Probe base incorporating indicator light
2. Eye recognition transmitter
3. Eye recognition sensor
4. Cheek support
5. Forehead support

TOP PARTS (Figure 2)

6. Measurement button
7. Forehead support position indicator

SIDE PARTS (Figure 3)

8. Forehead support adjustment wheel
9. Cheek support adjustment wheel

BOTTOM PARTS (Figure 4)

10. Battery cover
11. Silicone lid (USB cover)
12. Type label

BACK PANEL (Figure 5)

13. Measurement button
14. LOAD light
15. MEASURE light
16. REPEAT light
17. DONE light
18. SERVICE light
19. BATTERY light
20. POWER button

10. SETTING UP THE TONOMETER

Setting up your Icare HOME tonometer is easy, with few steps. The following subchapters describe how you get started.

10.1 INSTALLING OR CHANGING THE BATTERY

Lift the silicone lid that protects the USB port and keeps the battery compartment cover in place. Open the battery compartment cover by pressing the silicone lid slightly and sliding the battery compartment cover as shown in the Figure 6.

1. Silicone lid
2. Battery cover

Insert two CR123A lithium batteries in the correct order: (+) end upwards as shown in the figure left. Close the cover firmly and press the silicone lid in place to cover the USB port.

10.2 TURNING THE TONOMETER ON

Press the power button (20) to turn the tonometer on. The lights (14-19) are turned on briefly. Following a brief pause, the Load light flashes on the back panel to remind the user to load the single-use probe into the tonometer prior to measurement.

10.3 LOADING THE PROBE

The Icare HOME tonometer uses single-use probes that are packed in a plastic tube and wrapped in blister packs as shown in the Figures 8 and 9.

To load the probe:

1. Unwrap the probe (Figure 10).
2. Remove the lid of the probe container as shown in the Figure 11. Point the tonometer upward.
3. Drop the probe into the probe base by turning the probe container upside down as shown in the Figure 12.
4. Press the measurement button briefly (1 second) to activate the probe (Figure 13).
5. The probe moves rapidly back and forth (Figure 14).
6. If the green Measure light flashes, the probe is loaded correctly and ready for measurement (Figure 15).

11. USING THE TONOMETER**11.1 CHOOSING THE MEASUREMENT MODE**

The tonometer can operate in two modes:

Series mode

In the series mode, keeping the button pressed down (see Figure 16) initiates the measurement function, and the tonometer takes six rapid measurements

one after the other to obtain the final IOP reading. The button must be pressed for at least 3 seconds.

Single mode

You can use the single mode to take individual measurements one at a time. The single mode is especially useful for those patients who tend to blink heavily. Here you press the measurement button briefly (1 second) for each of the six measurements to obtain the final IOP reading (Figure 17).

11.2 AUTOMATIC EYE RECOGNITION

The tonometer includes an automatic eye recognition system that identifies which eye, right or left, you are measuring. The system has two infrared LED transmitters just below the probe base and one infrared LED sensor above the probe base, as in Figure 21. The right-hand transmitter sends invisible infrared light to the right and the left-hand transmitter to the left. The infrared light reflects from your nose to the sensor. The sensor knows from which transmitter the reflected infrared light came, and thus which eye you are measuring. The resulting eye indication is included in the data that you can transfer, as described in the section 11. Do not cover the eye recognition sensors, because covering the sensors causes an error.

EYE RECOGNITION COMPONENTS (Figure 21)

1. Left and right infrared transmitter
2. Infrared sensor

11.3 ADJUSTING THE MEASUREMENT POSITION

The tonometer has two adjustable supports (see points 4 and 5 in Figure 1), one for the forehead and one for the cheek, as shown in Figure 18. The supports are for ensuring accurate measurement distance and alignment.

To adjust the measurement position:

1. Adjust the supports using the adjustment wheels as shown in Figure 18.
2. Keep the probe horizontal and pointing perpendicularly to the center of the cornea.
3. Set the distance between the tip of the probe and the center of the cornea to be 4-8 mm (5/32-5/16") as shown in Figure 19.
4. Read the distance setting (forehead A•1, A•2, etc., cheek B•1, B•2, etc.) between the arrows on the scale (see point 7 in Figure 2) of the supports (see Figure 20) and write it down on a support position tag for a reference if the setting is changed. The support position tag is found in the Icare tonometer carrying case.
5. Do the same for the other eye unless only one eye needs monitoring.
6. Regularly verify that the support positions are correct. If the distances are incorrect during measurement, there will be error indications (section 13).

11.4 TAKING THE MEASUREMENTS

The probe will make a gentle and brief contact with the eye when you take the measurement. No topical anesthetic is needed. The recommended frequency of measurements is 3-4 daily with a maximum of 5-6.

To measure intraocular pressure:

1. Check that the MEASURE light still flashes on the back panel.
2. If the MEASURE light does not flash, press the power button and wait until the MEASURE light illuminates again.
3. Look straight ahead at a specific point while keeping eyes wide open as shown in Figure 22.
4. Bring the tonometer near the eye, the probe pointing perpendicular to the center of the cornea without a vertical or horizontal tilt. The position is correct when the probe base light is green and appears symmetrically in the center of your view. See Figure 23 and 24.
 - Correct alignment of the tonometer (see Figure 25).
 - Incorrect alignment of the tonometer (see Figure 26).
 - Incorrect alignment of the tonometer. Readjust so that you see only the front of the tonometer and the green light symmetrically in the center of your view (see Figure 27).

5. Press the measurement button:

Single mode:

Press the button briefly (1 second) and you hear a short beep, repeat it to take one measurement at a time till you hear a long beep and see the Done light (17) illuminated on the back panel (Figure 5).

Series mode:

Keep the measurement button down to obtain the sequence of six measurements till you hear a long beep and see the Done light illuminated on the back panel. The button must be pressed for at least 3 seconds (Figure 29).

6. If an error occurs, press the Measure button briefly (1 second) and continue the measurement. See also section 13 Troubleshooting.
7. If both eyes are measured repeat above steps under 11.3 and 11.4 using your other eye.
8. Press the power button for three seconds to turn the tonometer off. The probe slides out of the probe base. Note that the tonometer will power off automatically if unused for 3 minutes.
9. Remove and discard the used probe properly.

12. TRANSFERRING THE MEASUREMENT DATA

The tonometer stores information on every complete measurement sequence of six measurements. The stored information includes the calculated final eye pressure reading in mmHg, time and date of the measurement, identification of the eye (right or left) and the quality level of the measurement. Uploading is easy:

1. Start Icare software in your mobile device or PC (Figure 30).
2. Connect the tonometer to the mobile device or PC using the USB cable. The Load and Measure lights will flash. If no lights flash or the SERVICE and BATTERY lights flash, reconnect the USB cable.
3. The internal clock of the tonometer is automatically updated to the mobile device's or PC's time by the Icare software at this point (Figure 31).

More information about **Icare software**

<http://www.icaretonometer.com>

13. TROUBLESHOOTING

The tonometer automatically monitors and controls the measurement position and speed of the probe during the measurements, and indicates errors with sounds and lights. The following table instructs you in error situations and explains what the different lights and sounds mean. All the back panel indicator lights are presented in the list after the table.

Error light	Error sound	Reason	Action
Battery (Figure 32)	No	Battery is soon empty.	Prepare to change batteries.
Battery light is flashing	No	Battery is empty.	Change batteries.
Probe base light is solid red (Figure 33)	No	Too much vertical tilt.	1. Press the measurement button again to clear the error message. 2. Position the tonometer horizontally so that the probe base light is green.

Error light	Error sound	Reason	Action
Probe base light is flashing red and Measure light turns off (Figure 34)	Two long beeps	Tonometer is not correctly aligned or probe hit eyelid or eyelash.	<ol style="list-style-type: none"> 1. Press the measurement button again to clear the error message. 2. Set the probe perpendicular to the center of the cornea. 3. Adjust the forehead/cheek supports for correct distance. 4. Repeat measurement.
	Three long beeps	Tonometer is too far from eye.	
	Five short beeps	Tonometer is too close to eye.	
Repeat light is flashing and Probe base light is flashing red (Figure 35)	Two long beeps	<ol style="list-style-type: none"> a) Too much IOP deviation during the measurement, because the user did not keep the tonometer stable. b) Eye was not recognized. 	<ol style="list-style-type: none"> 1. Press the measurement button again to clear the error message. 2. Repeat the measurement. 3. Do not move the tonometer during the measurements, remove your hand or fingers from the infrared transmitters and sensor, move the patient's hair away from his/her temple side of the eye.
Service light is flashing and Probe base is flashing red (Figure 36)	Two long beeps	Damaged or dirty probe or probe base.	<ol style="list-style-type: none"> 1. Change the probe and/or probe base. 2. Contact provider of the tonometer or a certified service center.

The indicator lights on the back panel:

OFF (Figure 37)

1. Lights off

INITIALIZE (Figure 38)

2. LOAD (green icon)
3. MEASURE (green icon)
4. REPEAT (yellow icon)
5. DONE (green icon)
6. SERVICE (red icon)
7. BATTERY (red icon)

LOAD PROBE (Figure 39)

8. LOAD (green icon)

MEASURE (Figure 40)

9. MEASURE (green icon)

REPEAT (Figure 41)

10. REPEAT (yellow icon)

DONE (Figure 42)

11. DONE (green icon)

SERVICE (Figure 43)

12. SERVICE (red icon)

LOW BATTERY (Figure 44)

13. BATTERY (red icon)

EMPTY BATTERY (Figure 45)

14. BATTERY (red icon)

14. REPLACING THE PROBE BASE

Conditions during the use or transport of the tonometer can include lint, dust, heat, moisture, hair, grease or other debris that can get to the probe base and impair the probe movement. To resolve these conditions replace the probe base at least every 6 months. Change the probe, replace the probe base if the SERVICE light flashes.

Instructions for replacing the probe base:

- Turn off the tonometer
- Unscrew the probe base collar and put it in a safe place
- Remove the probe base by tilting the tonometer downwards and use your fingers to pull the probe base out of the tonometer
- Insert a new probe base into the tonometer
- Screw the collar in, to lock the probe base

15. CLEANING AND DISINFECTION

The forehead and cheek supports must be cleaned for each new patient. The health care professionals should use a wipe dampened with a 70-100% isopropyl alcohol solution. The patients should use a cloth or paper towel dampened with water. Do not immerse the tonometer in water or other liquid. The tonometer must not be immersed or cleaned using too much water.

16. ACCESSORIES, DETACHABLE PARTS AND OTHER SUPPLIES

Part number	Product Description	Weight	Dimensions
540	Probe base	4g	38 mm × 7 mm
TA022-118	TA022 narrow collar	2g	20 mm × 13 mm
560	Wrist strap	3g	270 mm × 10 mm × 10 mm
TA022-044	Carrying case	210g	270 mm × 135 mm × 60 mm
7179	Battery cover	3g	26 mm × 23 mm × 7 mm
571	Battery 3 V, CR123A	17g	35 mm × 17 mm
TA022-037	Support position tags	40g	70 mm × 41 mm × 13 mm
575	USB cable	23g	1 m
113	Probe Icare TP022, 50 pcs/box	55,14g	195 mm × 82 mm × 35 mm

17. LIFETIME

The expected service life of the device is 5 years. The maintenance described in the chapters 12-15 is required during the the expected service life. The shelf life of the probes in their intact original packaging is 3 years. Please check the shelf life (expiry date) from the probe label. We recommend that you inspect the device for mechanical and functional damage and the safety labels for legibility annually/every 12 months. The life time of the batteries is about 5000 measurements.

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18. TECHNICAL AND PERFORMANCE DATA

Type: TA022

Dimensions: approximately 11 cm x 8 cm x 3 cm

Weight: approximately 150g

Power supply: 2 x CR123 non-rechargeable batteries (make sure you use batteries with built-in PTC protection, for example Energizer Lithium Photo 123 3V CR123A)

Measurement range: 5-50 mmHg

Accuracy (95% tolerance interval relatively to manometry): $\pm 1,2$ mmHg (≤ 20 mmHg) and ± 2.2 mmHg (> 20 mmHg)

Repeatability (coefficient of variation): $< 8\%$

The serial number is located on the inside of the battery compartment cover. The lot number of the probes is on the side of the probe box and the blister packing.

There are no electrical connections from the tonometer to the patient. All the parts of the tonometer are applied parts and the tonometer has BF-type electric shock protection.

Operation environment:

Temperature: +10 °C to +35 °C
 Relative humidity: 30% to 90%
 Atmospheric pressure: 800 hPa - 1060 hPa

Storage environment:

Temperature: -10 °C to +55 °C
 Relative humidity: 10% to 95%
 Atmospheric pressure: 700 hPa - 1060 hPa
 It is recommended to cool the tonometer for an hour from the maximum storage temperature or to warm the tonometer for an hour from the minimum storage temperature before use.

Transport environment:

Temperature: -40 °C to +70 °C
 Relative humidity: 10% to 95%
 Atmospheric pressure: 500 hPa - 1060 hPa

Environmental restrictions for professional use include:

- Medivac vehicles or similar where vibration or noise levels are so high that the user cannot hear error signals.

Environmental restrictions for lay operators (patients):

- Environments where noise is so high that the user cannot hear the error signals.

Mode of operation: continuous

19. SYMBOLS

Caution

See operating instructions for more information

BF-type device

Single-use disposable

Serial number

Use by <date>

Manufacturer

IP22

Protected against insertion of fingers and will not be damaged or become unsafe during a specified test in which it is exposed to vertically or nearly vertically dripping water.

Class 1 LED product

This product meets the power requirements for a Class 1 LED product to IEC/EN 60825-1 (2001) under normal operating conditions and those of single fault failure.

Keep dry

Manufacturing date

Lot number


Sterilized using irradiation


Stand by

Do not discard this product with other household-type waste. Send to appropriate facility for recovery and recycling. EU WEEE (European Union Directive for Waste of Electronic and Electrical Equipment)

Storage environment

Transport environment

Temperature limits

Humidity limits

Atmospheric pressure limits

20. ELECTROMAGNETIC DECLARATION

Icare HOME is class B equipment and needs special precautions regarding EMC and needs to be installed and put into service according to EMC information provided below.

Guidance and manufacturer's declaration–Electromagnetic emissions		
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.		
RF emissions CISPR 11	Group 1	Icare HOME (TA022) is battery operated and uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby equipment.
RF emissions CISPR 11	Class B	Icare HOME (TA022) is suitable for use in all establishments, including domestic establishments and those directly connected to public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	NOT APPLICABLE	Power level of Icare HOME (TA022) is below standard requirement, batteries are not rechargeable.
Voltage fluctuations flickering emissions IEC 61000-3-3	NOT APPLICABLE	Icare HOME (TA022) batteries are not rechargeable.

Guidance and manufacturer's declaration–Electromagnetic immunity			
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transients/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable.
Surge IEC 61000-4-5	±1 kV for line(s) to line(s) ±2 kV for line(s) to earth	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable.
Voltage dips, short interruption and voltage variations on power supply lines IEC 61000-4-11	0 % UT for 0.5 cycle (1 phase) 0 % UT for 1 cycle 70 % UT for 25/30 cycles (50/60 Hz) 0 % UT for 250/300 cycles (50/60 Hz)	NOT APPLICABLE	Icare HOME (TA022) tonometer is not operational when connected to an external computer, Icare HOME (TA022) batteries are not rechargeable.

<p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>30 A/m</p>	<p>30 A/m</p>	<p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p> <p> WARNING: Sources of power frequency magnetic field should be used no closer than 15 cm (6 inches) to any part of Icare HOME (TA022), including cables specified by the manufacturer. Otherwise, degradation of the performance could result.</p> <p>The measurement method of the Icare HOME tonometer is based on magnetic induction and therefore an external magnetic field in line with the probe may prevent the measurement. In such case the tonometer will continuously ask to repeat the measurement. Situation can be solved either by removing the source of interference from the vicinity of the device or by performing the measurement in different location with no such interference.</p>
--	---------------	---------------	---

Guidance and manufacturer's declaration - Electromagnetic immunity			
Icare HOME (TA022) is intended for use and should be used in the electromagnetic environment specified below.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-Guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz - 80 MHz	3 V	 WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Icare HOME (TA022), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Interference may occur in the vicinity of equipment marked with the following symbol: 
	6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	6 V	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz - 2,7 GHz 80% AM at 1 kHz	10 V/m	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	380-390 MHz 27 V/m; PM 50%; 18 Hz	27 V/m	
	430-470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	28 V/m	
	704-787 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	
	800-960 MHz 28 V/m; PM 50%; 18 Hz	28 V/m	
	1700-1990 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	2400-2570 MHz 28 V/m; PM 50%; 217 Hz	28 V/m	
	5100-5800 MHz 9 V/m; PM 50%; 217 Hz	9 V/m	

21. TRAINING

All health care professionals, patients and caregivers should practice the use of the tonometer with this instruction manual or the patient guide and a training video, or receive training from a trained professional prior to performing tonometry.

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